



**Be Right™**



## Test strips free chlorine, 0-600 mg/L, 100 tests

**Nr. produs:** 2890200

**RON Preț (fără TVA):** Contact

Disponibil

Hach test strips are manufactured with a unique reagent pad at the tip, using proven chemistries based on standard reference methods. Just dip the strip, wait until color development occurs, and compare the color of the reacted pad to the chart printed on the bottle.

**Ideal for anyone performing semi-quantitative pre-screenings, spot-checks, or field tests**

**Simplify testing and get results in less than 1 minute, in the field or lab**

**Eliminate chemical handling and clean-up**

**Reduce guesswork by screening samples**

**\*Not for use in sea or brackish water**

---

### Specificații

Cei mai mici pași de incrementare:	Steps: 0, 25, 50, 100, 200, 400, 600
Domeniu de măsurare:	0 - 600 mg/L $\text{Cl}_2$
Nota de subsol:	*mg/L unless otherwise noted; ppb = $\mu\text{g/L}$ ; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Parametru:	Free Chlorine, High Range
Platforma:	Test Strip
Unitate:	100 tests