



Sampling bag, sterile with dechlorinating agent, 100 mL, 100 pcs

Nr. produs:	2075333
RON Preț (fără TVA):	Contact
Expediere în 2 săptămâni	

Sterile (certified), disposable, transparent polyethylene bags with trusted "Whirl-Pak" closure. Contains non-nutritive pill with 10mg sodium thiosulfate (will dechlorinate up to 10 mg Cl_2 in 100 mL sample) four 100mL fill-lines.

Specificații

Cantitate:	100
Capacitate:	177 mL ml Fill Line