



Test kit, phosphorus, orthophosphate (reactive) , model PO-19A

Nr. produs: 224801
RON Preț (fără TVA): Contact
 Nu este afișată data expirării

Produse periculoase



Simple, low-cost color disc method is the most accurate visual test kit method

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specificații

Domeniu de măsurare:	Triple Range
Domeniu de masurare :	0 - 40 mg/L PO ₄
Interval de măsurare standard:	0 - 0.8 mg/L PO ₄ ; 0 - 4 mg/L PO ₄ ; 0 - 40 mg/L PO ₄
Model:	PO-19A
Model test kit:	PO-19A
	(Includes filtration for turbid samples.)
Numar teste:	100
Nume metoda:	Color disc/Ascorbic acid
Parameter:	Phosphorus
Parametru:	Phosphorus, Orthophosphate (reactive) - As PO ₄ ³⁻
Platforma:	Disc culoare
Smallest increment:	0.02
Tipul carcasei:	F