



Test kit, phosphate, total ortho-/meta-, model PO-24

Nr. produs: 225001
RON Preț (fără TVA): Contact
 Nu este afișată data expirării

Produse periculoase



Kit provides measurement of ortho- (reactive) and meta- (condensed) phosphate (requires hot acid digestion).

Kit provides rapid measurement of ortho- (reactive) and meta- (condensed) phosphate (requires hot acid digestion).

Specificații

Domeniu de măsurare:	Triple Range
Domeniu de masurare :	0 - 40 mg/L PO ₄
Interval de măsurare standard:	0 - 0.8 mg/L PO ₄ ; 0 - 4 mg/L PO ₄ ; 0 - 40 mg/L PO ₄
Model:	PO-24
Model test kit:	PO-24
Numar teste:	100 Clear
	100 Turbid
Nume metoda:	Color disc/Ascorbic acid
Parameter:	Phosphate
Parametru:	Phosphate, total ortho-/meta- - As PO ₄ ³⁻
Platforma:	Disc culoare
Smallest increment:	0.02
	0.1
	1
Tipul carcasei:	F