



Hardness (Total & Calcium) Test Kit, Model HAC-DT

Nr. produs:

2063900

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

Produse periculoase



Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10 - 4000 mg/L. Contains approximately 100 of each (Total and Calcium) test.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Specificații

Conținutul cutiei:	Includes 6 reagents, digital titrator, delivery tubes, graduated cylinder, flask, manual and carrying case.
Domeniu de măsurare:	10 - 4,000 mg/L CaCO_3
Model test kit:	HAC-DT
Numar teste:	100 Ca
	+100 T
Nume metoda:	Digital Titrator/EDTA
Parametru:	Hardness, total & Calcium - As CaCO_3
Platforma:	Titrator digital
Smallest increment:	0.1-10
Tipul carcasei:	F

Conținutul cutiei

Includes 6 reagents, digital titrator, delivery tubes, graduated cylinder, flask, manual and carrying case.