



Dissolved oxygen test kit, model OX-2P

Nr. produs:

146900

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

Produse periculoase



Low dissolved oxygen levels in surface water systems usually indicate serious pollution. Adequate amounts are crucial for fish life, but conversely, DO must be excluded from boiler feed water to prevent corrosion. This test kit utilizes drop count titration and contains 100 tests.

Method / Chemistry: Drop count titration

Replacement reagent set PN 2439802

Utilizes drop count titration - a simple, economical titration method

Specificații

Domeniu de măsurare:	0.2 - 4.0 mg/L & 1 - 20 mg/L O ₂
Model test kit:	OX-2P
Numar teste:	100
Nume metoda:	Drop count titration/ Modified Winkler
Parametru:	Dissolved Oxygen - As DO
Platforma:	Drop Count Titration
Smallest increment:	0.2
Tipul carcasei:	D