



Test kit, sulfite, model SU-DT

Nr. produs: 2063300
RON Preț (fără TVA): Contact
 Nu este afișată data expirării

Produse periculoase



Using the Digital Titrator, this kit contains 100 tests with a range of 10-800 mg/L.

Digital Titrator and standardized titration cartridges provide laboratory accuracy for field titrations

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specificații

Cei mai mici pași de incrementare:	0.4-4.0**
Domeniu de măsurare:	4 - 800 mg/L SO_3^{2-}
Model test kit:	SU-DT
Numar teste:	100
Nume metoda:	Digital Titrator/Iodometric
Parametru:	Sulfite - As SO_3^{2-}
Platforma:	Titratore digital
Tipul carcasei:	F