



Test kit, alkalinity, model AL-DT

Nr. produs:

2063700

RON Preț (fără TVA):

Contact

Nu este afișată data expedierii

Produse periculoase



Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Determine phenolphthalein and total alkalinity by visual endpoint titration. Using the Digital Titrator, this contains 100 tests with a range of 10-4,000 mg/L.

Digital Titrator and standardized titration cartridges provide laboratory accuracy for field titrations

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specificații

Cei mai mici pași de incrementare:	0.1-10
Domeniu de măsurare:	10 - 4,000 mg/L CaCO ₃
Model test kit:	AL-DT
Nota de subsol:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Numar teste:	100
Nume metoda:	Digital Titrator/Sulfuric acid
Parametru:	Alkalinity - (P) & Total (MO)
	as CaCO ₃
Platforma:	Titratör digital
Tipul carcasei:	F