



Beaker, Polypropylene, Low Form, 250 mL, 6/pk

Nr. produs: 108076

RON Preț (fără TVA): Contact

Expediere în 2 săptămâni

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specificații

Capacitate:	250 mL
Diametru exterior :	83.6 mm with spout
Material:	Polypropylene, Low Form
Material corp:	Polypropylene
Unitate:	6 /pk