



Flask, volumetric, polypropylene, 100 mL

Nr. produs: 1406042

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Polypropylene flask and screw cap. Durable with excellent chemical resistance, with the same accuracy as glass. Permanent color graduation line with marking spot. Calibrated “to contain/deliver” at 20°C. The 200/203 capacity flask has dual marks useful in DO determination via Winkler Titration.

Specificații

Capacitate:	100 mL
Certificări:	Not Applicable
Material corp:	Polypropylene
Unitate:	1 each