



Be Right™



Hydrazine color disc test kit, model HY-2

Nr. produs: 184900

RON Preț (fără TVA): Contact

Nu este afișată data expirării

Produse periculoase



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Utilizing the color disc method with a range of 0.02-1 mg/L, this kit contains 300 tests.

Measures Hydrazine, a common oxygen scavenger in boiler water, using the simple color disc method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specificații

Domeniu de măsurare:	0 - 1 mg/L N_2H_4
Interval de măsurare standard:	0 - 1 mg/L N_2H_4
Model:	HY-2
Model test kit:	HY-2
Numar teste:	300
Nume metoda:	Color disc/ p-Dimethylamino- benzaldehyde
Parameter:	Detergents
Parametru:	Hydrazine - As N_2H_4
Platforma:	Disc culoare
Smallest increment:	0
Tipul carcasei:	D