



Tubing, tygon, R-3603, 4.8 mm (I.D.)

Nr. produs: 2074043

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Chemical-resistant, Formula R-3603 Tygon tubing with flex strength greater than rubber. Unaffected by oxidizing agents and compatible with virtually any laboratory chemical. Bends readily to sharp radius curves and will maintain strong grip on glass or metal. Highly polished bore prevents interior build-up while facilitating cleaning.

Specificații

Diametru exterior :	5/16 in
Formula (Text):	R-3603
Grosime conducte de perete:	1/16 in
I. D.:	3/16 in
Lungime:	50 ft
Material corp:	Tygon