



**Be Right™**



## Test kit, chloride, model CD-51

Nr. produs:

208601

RON Preț (fără TVA):

Contact

Nu este afișată data expedierii

Produse periculoase



Found in nearly all natural waters, chloride affects human taste above 250 mg/L. High levels inhibit plant growth, and many industrial processes need to limit chloride concentration. The model CD-DT, multi-range kit, is supplied with the demineralizer bottle for high-range sample dilution. Contains 100 tests.

**Utilizes drop count titration - a simple, economical titration method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Specificații

Domeniu de măsurare:	500 - 10000 mg/L & 5000 - 100000 mg/L Cl <sup>-</sup>
Model test kit:	CD-51
Nota de subsol:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Numar teste:	100
Nume metoda:	Drop count titration/
	Silver nitrate
Parametru:	Chloride - As Cl <sup>-</sup>
Platforma:	Drop Count Titration
Smallest increment:	500
	5000
Tipul carcasei:	B