



Bottle, BOD, black, 300 mL

Nr. produs: 2086047
RON Preț (fără TVA): Contact
Nu este afișată data expedierii

The standard for BOD analysis. Borosilicate glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stopper included (60 and 300 mL sizes). Serialized bottles have silk-screened, sequential

Specificații

Body Material:	Glass
Cantitate:	1 each
Capacitate:	300 mL
Material:	Borosilicate Glass
Package size UOM:	1pk