



**Be Right™**



## **Flask, Erlenmeyer, glass w/screw cap, 125 mL**

**Nr. produs:** 2089643

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

---

### **Specificații**

Capacitate:	125 mL
Material:	Glass Screw Cap
Unitate:	1 each