



Centrifuge Tube, Conical Bottom, 10/pk

Nr. produs: 2115839

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Centrifuge tube with conical bottom. Polycarbonate with molded-in 0.1mL graduations. For low-speed centrifugation up to 6000 x g. Autoclavable. The 15 mL size is also used with Settlometer Kit (P/N 21157-00). Supplied in packages of 10, without closures. Autoclavable. Nalge.

Specificații

Cantitate:	10 /pk
Capacitate:	15 mL
Lungime:	120 mm