



Tube, culture, 16 x 100 mm, pk/6

Nr. produs: 2275806
RON Preț (fără TVA): Contact
Nu este afișată data expirării

Cap sold separately: product number 2241106 (6 pieces).

Cap sold separately: product number 2241106 (6 pieces).

Cap has a firmly bonded rubber liner, providing a sure and highly inert sealing surface.

Specificații

Cantitate:	6 pk
Capacitate:	Capacity to Neck 11 mL
Diametru:	16 mm
Lungime:	100 mm