



Be Right™



Bottle, Sampling, Polycarbonate, 500 mL, 24/pk

Nr. produs: 2324579

RON Preț (fără TVA): Contact

Nu este afișată data expirării

Transparent polycarbonate bottle with cap features wide mouth (3.8 cm/1.48" diameter) opening for easy filling and square design for maximum space efficiency. Autoclavable with molded-in graduations (do not engage threads for autoclave cycle.)

Specificații

Body Material:	Polycarbonate
Cantitate:	24 /pk
Capacitate:	500 mL