



## Test kit, molybdate, model MO-LR

Nr. produs: 2359300

RON Preț (fără TVA): Contact

Nu este afișată data expirării

Produse periculoase



Molybdates are gaining over chromates for use as corrosion inhibitors due to lower toxicity. The high cost of these treatments warrant careful control of their concentration. Kit reagents react with Molybdate to form a yellow color, which is measured with a color disc.

**Kit reagents react with Molybdate to form a yellow color, which is measured with an easy-to-use color disc method**

**Effectively monitors corrosion inhibitors**

---

### Specificații

Domeniu de măsurare:	Single Range
Interval de măsurare standard:	0 - 3.0 mg/L Mo
Model:	MO-LR
Model test kit:	MO-LR
Numar teste:	100
Nume metoda:	Color disc/ Ternary complex
Parameter:	Molybdate
Parametru:	Molybdate - As Molybdenum
Platforma:	Disc culoare
Smallest increment:	0.25
Tipul carcasei:	D