



Filter, acrodisc, pore size 0.2 µm

Nr. produs: 2415600

RON Preț (fără TVA): Contact

Nu este afișată data expirării

For filtration of dilute biological fluids. Sterilized by gamma irradiation to eliminate potential contamination by EtO residuals. Inlet/Outlet connections: 13/25 mm: Female Luer-Lok inlet, standard male slip luer outlet. Pore size: 0.2 µm. Diameter: 25mm.

Specificații

Diametru:	NA
Dimensiune por:	0.2 µm
Filter Type:	Filter, acrodisc
Material:	Polysulfone
Packing:	NA