



## Filter holder, easy pressure syringe, PALL

**Nr. produs:** 246800

RON Preț (fără TVA): Contact  
Disponibil

Pack of 6 filter holders. Accepts 25 mm filter discs for fast filtration of 1-10 mL of sample.

**Fits most standard syringes with Luer-type fittings**

**Made of POM plastic with type-316 stainless steel membrane support screen**

**Autoclavable**

**Seal plugs included**

**Single filter PN 246801.**

---

### Specificații

Diametru exterior : 25 mm

Unitate: PK 6