



Rack, Test Tube, Polyethylene, 13 mm (O.D.), 90 Holes

Nr. produs: 2497900

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Ninety-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5").

Specificații

Cantitate:	1 each
Diametru:	13 mm
Dispunere eprubete:	6 x 15
Orificii eprubeta:	90