



Test tube, rack, polyethylene, 16 mm (O.D.), 60 holes

Nr. produs: 2497901

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Sixty-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5"). Color may vary.

Specificații

Cantitate:	1 each
Diametru:	16 mm
Disponere eprubete:	5 x 12
Orificii eprubeta:	60