



Test tube, rack, polyethylene, 25 mm (O.D.), 24 holes

Nr. produs: 2497903

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Twenty to ninety-place rack features a special blend of polypropylene which sinks in water, assuring non-floating racks when immersed. Autoclavable. Stackable and autoclavable. Overall dimensions: 24.6 x 10.7 x 6.4 cm (9.7 x 4.2 x 2.5").

Specificații

Cantitate:	1 each
Diametru:	25 mm
Disponere eprubete:	3 x 8
Orificii eprubeta:	24