



pH Paper, 0 - 14 pH Range, 100 tests

Nr. produs: 2601300

RON Preț (fără TVA): Contact

Disponibil

Provides simple, fast, economical method of pH measurement.

Provides simple, fast, economical method of pH measurement.

Provides a simple, fast, economical method of pH measurement. Obtain accurate results in weakly buffered, turbid and lightly colored samples.

Test Strips are great for quick evaluations of pH levels in the lab or in the field

Easy to use, disposable, and inexpensive!

Specificații

Cei mai mici pasi de incrementare:	0, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14
Domeniu de măsurare:	0 - 14 pH units
Numar metoda:	10056
Numar teste:	100
Parametru:	pH
Platforma:	Test Strip
Unitate:	100 tests