



## Cylinder, Graduated, 25 mL, Certified

**Nr. produs:** 2636240

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Certified Class A, this glass cylinder is calibrated "to deliver" according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in white marking. Bumper guard included for all sizes  $\geq 25$  mL.

---

### Specificații

Cantitate:	1 each
Capacitate:	25 mL
Certificări:	Class A
Înălțime:	192 mm
Material corp:	Glass
Metoda de măsurare:	To Deliver
Tip:	Simplu Metric