



Cylinder, Graduated, 500 mL, Certified

Nr. produs: 2636249
RON Preț (fără TVA): Contact
Nu este afișată data expedierii

Certified Class A, this glass cylinder is calibrated "to deliver" according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in white marking. Bumper guard included for all sizes ≥ 25 mL.

Specificații

Cantitate:	1 each
Capacitate:	500 mL
Certificări:	Class A
Înălțime:	390 mm
Material corp:	Glass
Metoda de măsurare:	To Deliver
Tip:	Simplu Metric