



Be Right™



Gas washing bottle, 1200 mL

Nr. produs: 2662200

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

For use in precipitating elements from gas. Used in the Hach Cold Vapor Mercury Apparatus Set (Prod. No. 2674400) to perform Cold Vapor Mercury Concentration Analysis—Hach Method 10065.

Specificații

Body Material:	Glass
Cantitate:	1 each
Capacitate:	1200 mL
Material:	Borosilicate Glass
Package size UOM:	1pk