



Iron (Ferrous) Color Disc Test Kit, Model IR-18C, 0-7 mg/L

Nr. produs:

2667200

RON Preț (fără TVA):

Contact

Expediere în 2 săptămâni

Produse periculoase



The most accurate visual test kit method, Hach color discs provide fast, accurate comparisons. Range: 0-7 mg/L. Contains 100 tests.

Simple, fast, color disc method

Rugged test kit contains all required reagents and apparatus to perform analysis in the field

Specificații

Domeniu de măsurare:	High Single Range
Model:	IR-18C
Model test kit:	IR-18C
Numar teste:	100
Nume metoda:	Color disc/ 1,10 Phenanthroline
Parameter:	Iron Ferrous
Parametru:	Iron, medium range - As Fe ²⁺ (ferrous iron)
Platforma:	Disc culoare
Smallest increment:	0.2
Tip:	Simplu Metric
Tipul carcasei:	B