



Thermometer, Non-Mercury, -30 - 120 °F

Nr. produs: 2676100

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Compact and durable, non-mercury thermometers enclosed in a protective nylon case. These thermometers are designed for ruggedness and come in many standard ranges. Suitable for all field and portable uses where an unprotected thermometer may break.

Specificații

Cantitate:	1 each
Diviziuni termometru:	2 °F
Domeniu de temperatură:	-30 °F - 120 °F
Tip fluid (termometru):	Non-Mercury