



## Thermometer, Non-Mercury, -5 - 45 °C

Nr. produs: 2676300

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Compact and durable, non-mercury thermometers enclosed in a protective nylon case. These thermometers are designed for ruggedness and come in many standard ranges. Suitable for all field and portable uses where an unprotected thermometer may break.

---

### Specificații

|                         |               |
|-------------------------|---------------|
| Cantitate:              | 1 each        |
| Diviziuni termometru:   | 0.5 °C        |
| Domeniu de temperatură: | -5 °C - 45 °C |
| Tip fluid (termometru): | Non-Mercury   |