



Thermometer, Dual Range, Non-Mercury, -10 - 110 °C, 0 - 220 °F

Nr. produs: 2676400

RON Preț (fără TVA): Contact
Nu este afișată data expedierii

Compact and durable, non-mercury thermometers enclosed in a protective nylon case. This thermometer is designed for ruggedness. Suitable for all field and portable uses where an unprotected thermometer may break.

Specificații

Cantitate:	1 each
Diviziuni termometru:	1 °C
Domeniu de temperatură:	-10 °C - 110 °C & 0 °F - 220 °F
Tip fluid (termometru):	Non-Mercury