



Test kit, hypochlorite, model CN-HRDT

Nr. produs:

2687100

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

Produse periculoase



Hypochlorite Kit, Model CN-HRDT with Digital Titrator

Hypochlorite Kit, Model CN-HRDT with Digital Titrator

Utilizing the Hach Digital Titrator, this kit is used to determine hypochlorite, high range (Bleach) solution concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 5 - 15% as Cl₂: 50 - 150 g/L as Cl₂. Contains 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Designed for on-site analysis of bulk hypochlorite deliveries or bulk storage solutions at drinking water or chloramination facilities

Verify concentration of hypochlorite solutions at time of delivery

Kit contains all required reagents, solutions, and glassware

Specificații

Domeniu de măsurare:	5 - 15 % = 50-150 g/L Cl ₂
Model test kit:	CN-HRDT
Numar metoda:	10100
Numar teste:	100
Nume metoda:	Digital Titrator/ Thiosulfate
Parametru:	Hypochlorite, high range (Bleach) - As Cl ₂
Platforma:	Titrator digital
Smallest increment:	0.05%
Tipul carcasei:	G