



Bottle, storage, autoclavable, wide mouth, 125 mL pk/12

Nr. produs: 2939103

RON Preț (fără TVA): Contact

Nu este afișată data expirării

Polyethylene bottles suitable for sample collection to chemical storage. Matched polypropylene caps have a valve-seal design for a perfect, leakproof fit. A protruding ring under the cap holds shrink-wrapping securely in place. Autoclavable 121°C (Do not autoclave with caps in place). Available in Narrow Mouth or Wide Mouth.

Specificații

Body Material:	Plastic
Cantitate:	12 /pk
Capacitate:	125 mL
Material:	PP
Package size UOM:	12pk