



Beaker, griffin, low form, glass, 1000 mL

Nr. produs:	50053
RON Preț (fără TVA):	Contact
Disponibil	

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologic

Specificații

Capacitate:	1000 mL
Material:	Glass, Griffin, Low Form
Material corp:	Glass
Unitate:	1 each