



Beaker, griffin, low form, glass, 2000 mL

Nr. produs: 50054
RON Preț (fără TVA): Contact
Nu este afișată data expirării

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologic

Specificații

Capacitate:	2000 mL
Material:	Glass, Griffin, Low Form
Material corp:	Glass
Unitate:	1 each