



Beaker, griffin, low form, glass, 100 mL, 12/pk

Nr. produs: 50072H

RON Preț (fără TVA): Contact

Nu este afișată data expirării

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologic

Specificații

Capacitate:	100 mL
Material:	Glass, Griffin, Low Form
Material corp:	Glass
Unitate:	12 /pk