



## QL filters, pore size 8 - 12 $\mu\text{m}$ , diameter 90 mm, circular, 100/ pk, qualitative

**Nr. produs:** 50655  
**RON Preț (fără TVA):** Contact  
Nu este afișată data expirării

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

---

### Specificații

Cantitate:	100 /pk
Diametru:	90 mm
Material:	Cellulose