



**Be Right™**



## Funnel, separatory, glass, 500 mL

**Nr. produs:** 52049

RON Preț (fără TVA): Contact

Nu este afișată data expirării

Pear-shaped (Squibb) separatory funnel with PTFE stopcock. Clear glass body allows for easy detection of layers between immiscible liquids, all the way to the stopcock. Highly chemical resistant and fully autoclavable. High density linear polyethylene s

---

### Specificații

|             |        |
|-------------|--------|
| Cantitate:  | 1 each |
| Capacitate: | 500 mL |