



Glass Fiber Filter, Diameter 90 mm, Preweighed for Standard Methods TSS procedure, 100/pk

Nr. produs: 2546140

RON Preț (fără TVA): Contact

Nu este afișată data expirării

Individually preweighed filter papers, preweighed in compliance with Standard Methods Total Suspended Solids procedure. Each in a planchet with tare weight listed. Each tare weight meets Standard Methods specifications.

934-AH is a registered trademark of Whatman, Inc./GE Healthcare

Specificații

Cantitate:	100 /pk
Diametru:	90 mm
Dimensiune por:	1.5 μ m
Material:	Glass Microfiber
Unitate:	100/pk